Multi-parameter tester



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

- Measures pH, conductivity, TDS, salinity, and temperature
- Graphic LCD display
- Switches between °C and °F
- Replaceable electrode





+961 70 913 590

sales2@mefosa.com



www.mefosa.com

pH/EC80 (Liquid)



Advanced pH meter

- pH and temperature tester
- Temperature range 0 to 99.9°C
- Graphical LCD
- Replaceable electrode





+961 70 913 590

sales2@mefosa.com



pH/ 618N (liqiud)



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال اقريقيا

JENCO







pH/6219N (Liquid)

- Two LCD display pH, mV, direct ion concentration and temperature
- Temperature range: -5.0 to 125.0°C
- RS-232 computer interface with software
- AC adaptor or battery power (low battery indicator





+961 70 913 590

sales2@mefosa.com



WATER TESTING STRIPS



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

- Free chlorine & pH strips
- Chlorine dioxide strips
- Free chlorine strips
- Extended range pH strips



2			5	
Wede	hlorine	Piks" 1		
They free Child	whe from 0-25 p	pre vent	14	
And and a second	* Dest accuracy of	oraul the		
State of the local division of the local div	Similar onion down for calibring agent in any present in the Colin Loss and any Similar processing of the Similar procesin	a and a	12	
C	SERENTEZON	-)	



	SonSafe 2
1	REE CHLORINE & PH
	STOCKLOSED STOCKLOSED
	TIT PROCEDURE:
N	Dip one test strip into a 50ml (about 2oz) water sample to
2	10 seconds with a constant. pentle, back and forth motion.
į,	Femove the strip and shake
	once, briskly, to remove excess valer.
	Mait 10 seconds.
	Using the color chart, match pH and Free Chilorine in this
)	order. Complete color
'	matching within 30 seconds.



+961 70 913 590

sales2@mefosa.com





Pool and Spa Test Strips

Measures:

- Total Chlorine
- Bromine
- Total Hardness
- Free Chlorine
- Total Alkalinity
- pH
- Cyanuric Acid





+961 70 913 590

sales2@mefosa.com







Metal Detector

- Suitable for most of packed or unpacked food
- Rejection feedback, air pressure monitor, fill up detection and locked bin optional
- Records could be export to USB drive (CSV format)





+961 70 913 590

sales2@mefosa.com





Sterilized Sampling Bags




MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

Provide a secure, contaminant free flexible container that ensures dependable analysis results.

A proven, economical and efficient way to collect, store and transport samples.









+961 70 913 590

sales2@mefosa.com



Moisture Analyzer



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

- Provides fast and accurate results
- High accuracy
- User-friendly touch-button design for easy operation
- Halogen heating and dual weight-temperature calibration

SPECIAL DRI





+961 70 913 590

sales2@mefosa.com



Incubator



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

Air inlet and outlet on the back, the air volume can be adjusted.

Adopt silicone rubber sealing strip and double glass door observation window, easy and clear observation.

Microcomputer controller, precise and reliable temperature control.







+961 70 913 590

sales2@mefosa.com





Digital Viscosity Meter

- Portable and compact size
- Can be powered by battery
- Only small volume of sample is required to measure







+961 70 913 590

sales2@mefosa.com



Heat Seal Tester



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائية للشرق الأوسط وشمال افريقيا

Heat Seal Tester is composed of upper and lower heatsealing jaws.

Before the test, preset the heat seal temperature, pressure and dwell time value, place the specimen in between the upper and lower jaws, and then press start button.

The whole sealing process can be finished automatically.





+961 70 913 590

sales2@mefosa.com

REFRACTOMETER



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائية للشرق الأوسط وشمال افريقيا

Refractometers measure the refractive index of liquids, pastes and solids with high precision

- Digital handheld refractometer
- Handheld
 refractometer
- Digital refractometer
- Abbe refractometer
- Process refractometer







+961 70 913 590

sales2@mefosa.com





Temperature, pressure and humidity data loggers

Monitoring solutions for the food & beverages, and catering industry.

The systems allow accurate and specific control of the quality of conservation, refrigeration, heat treatments and cooking.





+961 70 913 590

sales2@mefosa.com







Digital Thermometer

Use the Dickson digital thermometer to ensure that the food is cooked to safe temperatures, killing harmful pathogens and preventing foodborne illnesses.



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائية للشرق الأوسط وشمال افريقيا









www.mefosa.com



sales4@mefosa.com



DISPLAY USB TEMPERAURE LOGGER

- Temperature range: 0-70°C
- It is easy to collect temperature and humidity data
- Replaceable sensor







+961 70 913 590

sales2@mefosa.com





TEMPERATURE Report Logger

- Temperature range: -34 to 80°C
- IP63 Case
- Simple, succinct graph
- Long battery life
- No software required







+961 70 913 590

sales2@mefosa.com





Pathogen Q-PCR Detection Kits

- Results can be obtained within 24 hours.
- Detect 1–10 CFU/25 g sample, including enrichment
- Includes DNA extraction and real-time PCR analysis
- Licensed





+961 70 913 590

sales2@mefosa.com





Rapid Test Kits

- Algal toxins
- Bacterial toxins
- Industrial chemicals
- Mycotoxins
- Glyphosate test kits
- Estrogens
- Veterinary residues and additives



- Pesticides and herbicides
- GMOs



()



+961 70 913 590

sales2@mefosa.com

Portable Water Quality Laboratory

Kit includes:

- 2 Wagtech incubators
- Digital Arsenic Test System
- Lumiso Expert
- Turbimeter Plus
- Micro 800 pH/Temperature Meter
- Micro 800
 Conductivity/TDS meter



+961 70 913 590

sales2@mefosa.com



SPECI

www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

Pool Check® 7-Way Test Strips

Measuring for 7 of the most important pool & spa parameters.

Suitable for pool and spa





+961 70 913 590

sales2@mefosa.com



www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا



Chlorine Dioxide High Test Strips

After chlorination treatment, measures critical Chlorine Dioxide in water ranges (from 10– 2500 ppm) - NO BLEACH OUT of indicator like other

tests.

Suitable for pool and spa





+961 70 913 590


sales2@mefosa.com



www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا



SenSafe® Free Chlorine Test Strips

Perfect for determining whether or not you have the correct amount of sanitizer present.

Suitable for pool and spa





+961 70 913 590

sales2@mefosa.com



www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا



SenSafe® Chlorine/pH

SenSafe® Chlorine 2 way is the ideal choice for measuring Free Chlorine and pH in water.

Suitable for pool and spa





+961 70 913 590

sales2@mefosa.com



www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا



Pool Check® 6-Way Test Strips

Each strip tests:

- total chlorine
- total hardness
- free chlorine
- total alkalinity
- pH
- bromine

Suitable for pool and spa



+961 70 913 590





www.mefosa.com



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا





eXact iDip® Photometer

- Handheld photometerWaterproof
- Built-in sample cell
 Free downloadable app with pre-installed tests: pH, free chlorine, total chlorine, and total alkalinity



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا





Suitable for pool and spa



+961 70 913 590





Lumitester Smart



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائية للشرق الأوسط وشمال افريقيا

Portable luminometer for ATP test. The ATP Test provides easy inspection for cleanliness levels by using LuciPac







+961 70 913 590

sales2@mefosa.com



Common Pool Water



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا

Issues

Algae growth
Chemical Imbalance
Contamination
Filtration Problems







Lumitester Smart



MENA FOOD SAFETY ASSOCIATES مستشاره السلامة الغذائية للشيرة الأوسط وشمال افريقيا

Portable luminometer for ATP test. The ATP Test provides easy inspection for cleanliness levels by using LuciPac

Suitable for pool and spa







+961 70 913 590

sales2@mefosa.com



Electronic Precision Balance



- Repeatability ± 0.001 g
- Response time: 3-5 seconds



MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائية للشرق الأوسط وشمال افريقيا







sales2@mefosa.com



Candy / Deep Fry Thermometer

- Temperature range: 40-200°C
- Large, easy grip insulated handle with built-in pan clip.
- Adjust immersion depth with the included clip





MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا





+961 70 913 590

sales2@mefosa.com





Electronic Balance

- 5100g capacity
- 0.1 g readability
- Soft backlight LCD display
- Built-in RS232 interface







+961 70 913 590

sales2@mefosa.com





Digital Refractometer

- Brix 0-45%
- Refractive index 1.3330 -1.4098
- Temperature range: 0-40°C









sales2@mefosa.com






Handheld Salinity Meter

- Measures the amount of sodium chloride and temperature
- TDS range: 1.0 to 50.0 ppt
- Temperature range: -10.0 to 110.0
 °C







sales2@mefosa.com





LabStart

- Inexpensive
- Requires a power supply
- üBread and bread products
- üCereals
- üPowder, spices, and dried products
- Accuracy: +/-0.030 aw within calibration range
- Calibration points: 33%, 75%
- Temperature range: 15-30°C



novasina



+961 70 913 590







LabSwift

- More advanced
- Can operate on battery power
- üBread and bread products
- üCereals
- üPowder, spices, and dried products
- üMilk powder, dry dairy products
- Accuracy: +/-0.010 aw within calibration range
- Temperature range: 5-45°C



novasina



+961 70 913 590



sales2@mefosa.com







VS



LabStart



+961 70 913 590

sales2@mefosa.com







COMPACT DRY

Compact Dry is a ready-touse test method which helps to reduce the time needed to perform microbial testing. Therefore, it allows maximum productivity by increasing efficiency. The plates can be used to test raw materials as well as finished products like food, beverage, meat, cosmetic or other samples.





+961 70 913 590





(∰)



Easy Plate

A convenient, prepared media plates for faster and easier microbiological testing.

Easy Plate comes in an easy-to-use, stackable film plate format that saves on space, reduces waste, eliminates preparation time and labor, and streamlines your test procedures





+961 70 913 590

sales2@mefosa.com




(曲)



Food Stamp

Prepared culture media for self-hygiene control at all kinds of food related industries.





sales2@mefosa.com





LISTERIA SWAB

Developed to reduce the risk of contamination on work surfaces or in finished products by Listeria spp., but especially by L. monocytogenes.

The test is based on the ability of Listeria spp. to hydrolyze esculin, in the medium enriched with esculin, this hydrolysis leads to the formation of an unmistakable, black-colored precipitate (precipitate).







+961 70 913 590

sales2@mefosa.com



COLIFORM SWAB

Easy to use early warning test for coliform contamination of work areas.

- Simple color change indicates contamination
- No color change mean's its clear
- Self-contained swab test. No other reagents to add
- Results can be available within 18 hours
- Can detect as low as 1cfu per cm^2
- Ideal for HACCP program.
- Production is ISO 9007 approved



+961 70 913 590





MENA FOOD SAFETY ASSOCIATES مستشارو السلامة الغذائبة للشرق الأوسط وشمال افريقيا



sales2@mefosa.com () ww



Labtouch-aw

A semi temperature-control as well as integrated surface temperature measurement allows an accurate and reproducible measurement.





sales2@mefosa.com





Labmaster neo-aw

Outstanding measuring instrument for effortless and worry-free water activity measurement and documentation, including:

- 21CFR11 compliant audit trail and user management
- automatic verification or adjustment
- contactless salt standard identification by RFID
- fast mode that guarantees results within 10 minutes
- full thermostatisation in the range of 0-60°C
- measuring range from 0.0300-1.000aw
- a simple protective filter handling by magnet fixing

Of course also conformable with AOAC, ISO, USP and other standards.





+961 70 913 590







ANALYTICAL BALANCES





sales2@mefosa.com





PRECISION BALANCES






sales2@mefosa.com



INSPIRATIONAL MOOD

Creating a fashion moodboard is a fun and creative way to explore your personal style and inspiration. Whether you're looking to revamp your wardrobe or just want to experiment with different looks, a moodboard can help you visualize your ideas and bring them to life.



Handheld Salinity Meter

- Measures the amount of sodium chloride and temperature
- TDS range: 1.0 to 50.0 ppt
- Temperature range: -10.0 to 110.0 °C



- +961 70 913 590
- sales2@mefosa.com (www.mefosa.com

Electronic Balance

- 5100g capacity
- 0.1 g readability
- Soft backlight LCD display
- Built-in RS232 interface







Digital Refractometer

Brix 0-45%

SPECIAL

- Refractive index 1.3330 -1.4098
- Temperature range: 0-40°C



- +961 70 913 590

sales2@mefosa.com (www.mefosa.com

Candy / Deep Fry Thermometer

- Temperature range: 40-200°C
- Large, easy grip insulated handle with built-in pan clip.
- Adjust immersion depth with the included clip





+961 70 913 590



